




Microworld®

MODEL 2010/M PRISM COUPLER

BULK MATERIAL OR THICK FILM INDEX/BIREFRINGENCE MEASUREMENT

The Model 2010/M can operate as a fully-automated refractometer, providing high accuracy measurement of refractive index and index anisotropy for solid or liquid bulk materials and thick films in the index range 1.0 - 2.6 without use of toxic or corrosive index matching fluids. This bulk material index measurement is not an option, but a standard feature provided with every system.

For bulk material index measurements, if a material with index n (see Fig. 1) is brought into intimate contact with a prism with index n_p .



Laser light directed onto the base of the prism will be totally reflected at the prism base until the angle of incidence θ becomes less than the critical angle θ_c where

$$\theta_c = \arcsin(n/n_p) \quad (1)$$


In the Model 2010/M, the prism, sample, and coupling head are all mounted on a rotary table which can be turned under computer control to change the angle of incidence θ to scan through the critical angle for the prism/sample interface. θ_c is easily measured using the detector in Fig. 1 since the detector intensity drops abruptly as θ drops below the critical angle and light starts to leak into the bulk material. Since n_p is known, n can easily be determined from Equation (1).

This same approach can also be used to measure refractive index for thick films. In the case of a film, the critical angle of the prism/film interface establishes an upper limit to the angle at which film propagation modes can occur. As a film with fixed index gets thicker, the angle at which the first propagation mode occurs asymptotically approaches the critical angle as defined by equation (1) where n now refers to the film index rather than to a bulk index. If we make the approximation that the angle of the first film mode equals the critical angle, the error in the measured index due to this approximation is less than .004 for a film thickness of 3 microns, less than .001 at 5 microns, and less than .0003 at 10 microns. Thus, this approach can be applied with good accuracy to measure the index of materials with thickness ranging from several microns on up to "bulk" thickness.

For measurement of liquid samples, a cell to bring the liquid in direct contact with the prism is available.

The bulk index approach can also be used to study index anisotropy in solid materials. In normal operation (Fig. 2a), the Model 2010/M uses a polarized laser with TE incidence, i.e., electric field vibrating transverse to the plane of incidence (out of the plane of Fig. 1). Thus, if the material in figure 1 is simply manually

rotated against the prism (Fig. 2b), the electric field vector can be made to assume any orientation within the plane defined by the surface of the material, and the index for electric field vibration along any axis lying parallel to the material surface can be measured.



In addition, with option #2010/M-TM, which allows TM incidence (magnetic field transverse to, and electric field parallel to, the plane of incidence), index for electric field vibration perpendicular to the material surface is measurable (Fig. 2c). This is because the electric field vibrates perpendicular to the material surface at the critical angle since the refracted ray must propagate parallel to the surface at the critical angle.

Accuracy and resolution: With the high resolution rotary table (a no-cost option), worst case index accuracy is ± 0.001 and resolution is ± 0.0001 for bulk materials and films thicker than 10 microns. Absolute accuracy can be improved to approximately ± 0.0001 , however, if the user is willing to perform a simple calibration procedure (measurement of an absolute index standard) with each prism. Accuracy decreases with film thickness, with measured index systematically low by 0.001 at 5 microns, and by 0.003-0.004 at 3 microns.

Measurement time: Typical measurement time is 10-20 seconds.

Index measuring range: Minimum and maximum index measurable varies with prism type selected: 1.55-2.45 (200-P-2 prism), 1.25-2.10 (200-P-3), 1.20-2.00 (200-P-4), and 1.00-1.80 (200-P-1). Special 200-P-2 prisms are available to extend measuring range up to approximately 2.65.

Dispersion (index vs wavelength): Unlike most conventional refractometers, which are single-wavelength (typically 589 nm), the 2010/M can be equipped with as many as five lasers, allowing easy measurement of dispersion across a wide wavelength range. After index is measured at three or more wavelengths, Model 2010/M fitting software will generate extremely accurate curves of index vs wavelength in only a few seconds.



GRENOBLE - FRANCE

Tel : +33 (0)4 76 56 16 17

Email : contact@microworld.eu

www.microworld.eu